

OPTIM[®] 33TB

ONE-STEP CLEANER
& INTERMEDIATE DISINFECTANT

Active Ingredient: **0.5% Hydrogen Peroxide**



The Perfect Balance

No compromise between safety and efficacy.

OPTIM 33TB is one of the fastest and most effective cleaner & intermediate level disinfectants. OPTIM 33TB is safe; it does not require any GHS hazard pictograms, signal words or statements. OPTIM 33TB provides you with confidence that you are compliant and protecting your patients, staff, equipment and the environment; the perfect balance.



Effective Cleaner

- › Removes organic debris while disinfecting, saving you time and money¹



Fast Broad-Spectrum Disinfectant

- › 1 minute tuberculocidal, bactericidal, virucidal contact time
- › 3 minute fungicidal contact time
- › Ensures compliance, saving you time and money²



Safe for yourself and your patients

- › Non-irritating to skin, eyes and respiratory tract
- › Non-toxic to reproduction; child development and fertility
- › No VOCs (volatile organic compounds)



Environmentally-Friendly

- › Active ingredient, Hydrogen Peroxide, breaks down into water and oxygen
- › No active residues, reducing environmental impact




Time and Cost Efficient

- › Half as many wipes as the leading competitor are needed to keep your surface wet long enough to achieve disinfection³



Compatible

- › Designed to be compatible with a wide range of materials⁴



Fast Contact Times

A contact time is the amount of time a surface must stay wet for disinfection to occur. OPTIM 33TB's contact times are faster which ensures compliance and provides confidence that disinfection has been achieved.

OPTIM 33TB has been proven effective against difficult to kill pathogens such as Norovirus, giving you confidence that your facility is protected against pathogens of concern.

Surface Disinfectant Contact Time*

Product Name	Active Ingredient	Mycobacteria		Non-enveloped Viruses					Enveloped Viruses				Bacteria				Fungicidal
		TB	Adeno	Polio	Noro	Rota	RSV	HIV	HBV (Hep.B)	HCV (Hep.C)	Influenza A	MRSA	VRE	E.coli	PSA		
OPTIM 33TB	Hydrogen Peroxide	1		1	1	1		1	1	1	1	1	1	1	1	3	
Birex SE	Phenolics	10	10					10				10		10	10	10	
Patterson pdCare	Quats-alcohol	3						2	2	2	2	2	2	3	3	3	
Micro-Kleen 3	Quats-alcohol	3		3	3	3		3	3	3	3	3	3	3	3	3	
Caviwipes	Quats-alcohol	3						3		3	3	3	3	3	3	3	
Sanicloth Super	Quats-alcohol	2	2			2	2	2	2	2	2	2	2	2	2	2	
Germicide 3	Quats-alcohol	3		3				3	3	3	3	3	3	3	3	3	
Biotext	Alcohol Anhydrous	5													3	3	
BioSurf	Alcohol Anhydrous	1													20 sec	10 sec	
Tb Minuteman	Alcohol Anhydrous	1		1				1			1	1	1	1	1	1	
Clorox HP Wipes	Hydrogen Peroxide	5			3	1	30 sec	30 sec	30 sec	30 sec	1	1	30 sec	30 sec	30 sec	5	

* Contact time in minutes. If blank, no published information was available. Information was obtained from publicly available material in October of 2016. Information will be updated if there are changes in claims.

Ordering Information for OPTIM 33TB

Ready-to-use Pre-saturated Wipes

SCI33W12 OPTIM 33TB RTU Disinfectant Wipes, Unscented (12 x 160 sheets)

SCI33W12-ALOE OPTIM 33TB RTU Disinfectant Wipes, Aloe Scent (12 x 160 sheets)

Liquid

SCI334X4 OPTIM 33TB RTU Disinfectant, Unscented, 4L Bottle (4x4L)

SCI334X4-ALOE OPTIM 33TB RTU Disinfectant, Aloe Scent, 4L Bottle (4x4L)

SCI3312X1 OPTIM 33TB RTU Disinfectant, Unscented, 1L Bottle (12x1L)

SCI3312X1-ALOE OPTIM 33TB RTU Disinfectant, Aloe Scent, 1L Bottle (12x1L)

Dry Wipes Refill & Accessories

SCI33R18 Dry Wipes Refill, 6" x 7" (18 x 160 sheets / case)

SCI33R18XL Dry Wipes Refill, Large, 10" x 10" (18 x 60 sheets / case)

FSA01 Sprayer for 1L Bottle



For more information contact:

ANDA Medical

330 Marwood Drive, Unit 3, Oshawa, ON L1H8B4

Phone: 905-725-ANDA(2632)

sales@andamedical.com

¹ John A. Molinari, Ph.D., "OPTIM 33TB Cleaning Study," Ann Arbor: Dental Advisor, 2010.
² John A. Molinari, "Environmental Surface Wetness Test: Comparison of Disinfectant Wipes," Ann Arbor: Dental Advisor, 2015.
³ John A. Molinari, "Environmental Surface Wetness Test: Comparison of Disinfectant Wipes," Ann Arbor: Dental Advisor, 2015.
⁴ John Molinari, Ph.D., Peri Nelson, B.S. Effect of OPTIM 33TB on Tensile Properties of Dental Chair Fabrics. Ann Arbor: Dental Advisor, 2016.

