

INFRARED FOREHEAD THERMOMETER



- 1 second measurement
- Measurement in Celsius or Fahrenheit
 - +/- 0.4°C (32.0°C – 33.9°C)
 - +/- 0.3°C (34.0°C – 43.0°C)
- Turn on or off buzzer/speaker/vibration function
- Orange LED lighting, high fever alarm
- Automatic Shutdown
- Requires 2 x AAA batteries (not included)

Product Code: KS-IRFT

Minimum Order: 40ea (1 Carton)

Price: \$90.00 CAD/each (\$3,600 CAD/ Carton)

Prices as of 7/24/2020. Subject to change without notice